



## A NEW HPLC-DAD METHOD FOR STEROLS ANALYSIS IN BIOLOGICAL FLUIDS

Amaral C. (1,2)\*, Rodrigues M. R. (1), Mota Freitas M. (1), Lopes E. (1), Tomaz C. (2), Cardoso M. L. (1)

(1) Centro de Genética Médica Dr. Jacinto Magalhães, INSA

(2) Centro de Investigação em Ciências da Saúde (CICS), Universidade da Beira Interior.

**Introduction:** Several methods were described in literature for sterols quantification. Among them gas chromatography-mass spectrometry (GC-MS) with selected-ion monitoring (SIM) is the most common one. We set up a new sensitive and specific reversed-phase high-performance liquid chromatographic method using DAD detection (HPLC-DAD) for the simultaneous quantification of sterols in biological fluids as well as cell homogenates.

**Material and Methods:** We used HPLC-DAD equipment (Lachrom, Merck Hitachi®) with software HSM-system 1. After an injection of 100 µl, sterols were separated at 30°C, in a nucleosil C-18 column (250mm × 5 mm; 5µm particle size) with an isocratic elution. Mobile phase was methanol /water (100:4, v/v) and detection was performed at a wavelength of 210 nm. The flow rate used was 1 ml/min. and total running time was 40 min..

**Results:** We were able to identified ergosterol, 7-dehydrocholesterol, brassicasterol and cholesterol at retention times 25.41 min., 27.99 min., 28.77 min. and 32.15 min., respectively. The method was linear within the concentration range, 7.81–250 µM ( $r^2 > 0.9996$ ). The lower limit of quantification is 1nM for cholesterol, 7-dehydrocholesterol, ergosterol and brassicasterol. Within-run and between-run coefficients of variation were 1.65% and 9.99% for ergosterol and 2.55% and 10.54% for brassicasterol, respectively.

**Conclusion:** Our method is suitable for identification, separation and quantification of ergosterol, 7-dehydrocholesterol, brassicasterol and cholesterol in biological fluids and in cell homogenates. It was noted that our HPLC method had a better performance (for low sterol concentration) than GC-MS. Thus, this method is simple and accurate enough to be used in routine sterols analysis, and represents an interesting alternative to other methods described in literature.